

DICOM  
configuration

How to use

Copyright file

DICOM  
information



# CARE Analytics – Enhanced Dose Reporting

Operator Manual

Answers for life.

**SIEMENS**

# DICOM Configuration of CARE Analytics

**Note:** CARE Analytics is freeware and will not be installed by Siemens service engineers. The selection of a PC or USB stick in the hospital network is the responsibility of the hospital IT administrator. He has to ensure that there is no risk related to the introduction of this PC or USB stick in the DICOM network (i.e. virus protection). It is prohibited to put CARE Analytics on any Siemens medical device.


## Step 1

**Configure CARE Analytics' AE title and IP port.**

**Explanation:**

Worked as a DICOM SCU (Service Class User), this AE title and IP port is used to identify CARE Analytics.

- The AE title should not exceed 16 characters.
- This configuration part should not be changed frequently because the information is used by the DICOM node server. Once these parameters are changed, all the DICOM node servers that are listed in the "DICOM node list" will also need to be changed accordingly.




## Step 2

**Input DICOM node server's information (CARE Analytics' node name, AE title, node port, and IP address).**

**Explanation:**


For detailed configuration, please refer to the DICOM node server's configuration manual.

The following example of how to configure the DICOM node is provided for demonstration purposes:  
Server: Somaris/5 syngo CT 2009E(VX70A)  
IP Address: 139.24.185.19  
Client: Microsoft Windows XP  
IP Address: 139.24.185.236



# DICOM Configuration for CT Scanners


CARE Analytics is a stand-alone tool that can be used on any office computer. It can query/retrieve DICOM Structured Dose Reports from a DICOM node server. The retrieved structured reports can be exported to Microsoft® Office Excel for further analysis. This tool only supports retrieval and analysis of dose reports that conform to DICOM standard 2008 "CT RADIATION DOSE SR ID TEMPLATES". The templates that comprise the CT Radiation Dose SR are interconnected as depicted below:



Before starting a query on a new DICOM node server, both the DICOM node server and the CARE Analytics tool need to be configured.

## Step 1

Open the Service UI on the *syngo* system, select "Configuration → System Options". Enable checkbox "DICOM networking".




## Step 2

Configure entity titles for “Study Transfer”.


### Explanation:

1. “VB40” is just used as the DICOM node’s identification. The user can define it at will.
2. syngo’s DICOM node port is 104 (syngo default value and cannot be changed).
3. Target DICOM nodes IP address is “139.24.185.19”.




## Step 3

Add client’s IP address to Host list.



## Step 4

Add DICOM configuration, and then click "Save".



## Step 5

Click "Save" and "Finish", a restart of the application SW is necessary to make changes valid!

# DICOM Configuration of Artis zee Systems

CARE Analytics is freeware, a non-medical software tool designed to query dose information from Structured Dose Report objects stored on configured DICOM nodes. The tool can be executed from any PC in the hospital network or from a USB stick approved by a hospital IT administrator. CARE Analytics was tested with Structured Dose Reports from Artis zee/zeego systems beginning with version VC14.

## Preparation

### Step 1

**Alternative 1:**  
Report directly from the Artis zee system:  
Ask the customer for the following  
parameters for the PC to be connected.

Host Name	
TCP/IP Address	
Logical Name	
AE Title	

**Alternative 2:**  
Report from PACS: Ask the customer for  
the following parameters.

Logical Name	
--------------	--

**Explanation:**

A distinction must be made here whether the customer would like to have the report sent directly from the Artis zee system or from PACS.

### Step 2

**Note:** This point should be worked through only if the report is to come directly from the Artis zee system.

Host Name	
TCP/IP Address	
AE Title	

**Provide the following information to the customer:**

1. Open the service software.
2. Select "TCP/IP LAN" under "Local host".
3. Enter the name displayed under "Computer name" into the table.
4. Page to the "IP address" and enter the IP address below.






- 
- The screenshot displays the 'Siemens Med Service Software' interface within a web browser. The top navigation bar includes tabs for 'Configuration', 'Event Log', 'Reports', 'Documentation', 'Home', and 'Help'. On the left, a sidebar lists several configuration categories: 'System Options', 'Local Host', 'Security', 'Service', 'DICOM', 'Import/Export', and 'Artis System'. Each category has a list of sub-items. The main area is titled 'Identification' and contains several sections for network configuration:
- Computer name:** A text field containing 'axisec01'.
  - Adapter:** A dropdown menu set to 'DICOM : Intel(R) PRO/1000 PT Dual Port Server Adapter'.
  - IP Address:** Includes a radio button for 'Obtain an IP address from DHCP server' (which is selected), a text field for 'Specify an IP address', and a numeric input grid for IP address and subnet mask. The IP address is currently blank, and the subnet mask is '255.255.255.0'.
  - WINS:** Includes text fields for 'Primary WINS server' and 'Secondary WINS server', both currently blank, and checkboxes for 'Enable DNS for Windows Resolution' (unchecked) and 'Enable LMHOSTS Lookup' (checked).
  - DNS:** Includes a text field for 'Domain' (blank) and a numeric input grid for 'DNS service search order' (all digits are blank).
- At the bottom of the window, there are 'ready', 'Finish', and 'Save' buttons.

## Step 3

### Check the configuration for Dose Structured Report at the Artis zee system.


**Note:** This point must be checked, regardless of whether the report is to come directly from the Artis zee system or from PACS.

1. Open the service software.
2. Click on "Configuration".
3. Select the "Image System" point under "Artis Components".
4. Select "Configure Exam Protocol Transfer".
5. Click on "Config".
6. Check whether "Neither Dose SR objects or Exam Protocol XA image objects" or "Exam Protocol XA image objects" is selected.
  - If "Neither Dose SR objects or Exam Protocol XA image objects" is selected, select "Dose SR objects".
  - If "Exam Protocol XA image objects" is selected, select the "Dose SR objects and Exam Protocol XA image objects" point.
7. Click on "OK".
8. Click on "Save".





9. Select "Completed Flag Handling".
10. Click on "Config".
11. Check to ensure the point "Quant Reports / Example Protocol" is set to "enabled". If not, select "enabled".
12. Click on "OK".
13. Click on "Save".
14. Click on "Exit".
15. Click on "Finish".




## Configuration when the report is to come directly from the Artis zee system

**Note:** This point should be worked through only if the report is sent from the Artis zee system. This configuration is necessary to send the "Dose SR" to an extra DICOM node.


### Step 1

Set the DICOM configuration of the Artis zee system.

1. Open the service software.
2. Click on "Configuration".
3. Under "DICOM", select the "Network Nodes" point.




4. "Host Name" and enter the IP address under "TCP/IP Address". This information comes from the customer.
5. Click on ">".
6. Enter the AE Title under "edit Name". This information comes from the customer.
7. Select the "Host Name" entered under "Host" on the previous page.
8. Enter the same name under "Edit AE Title" as was entered under "Edit Name".
9. Enter the value "5104" under "Post Number".
10. Select "Storage".
11. Click on "Add".
12. Click on "Save".
13. Click on "OK".
14. Click on "Next".
  - A configuration check is performed automatically.
  - "Auto config needed" is displayed in the list.
15. Click on "Auto Conf".
  - When this is done, "Auto config" is started automatically.
  - "Auto Conf finished" appears in the status line.
16. Click on "Exit".



## Step 2

Check the Name under "Network Nodes".

1. Select the "Image System" point under "Artis Components".



## Configuration when the report is to come from PACS

**Note:** The definition for the Auto Report is checked here.


## Step 1

Select "Options". Select "Configuration"  
→ The Artis zee Configuration Panel is opened.



## Step 2

Double-click on "Transfer".



## Step 3

Check whether the PACS name is already present under "Destination".


If yes,

- Do not make any changes.
- Close the windows.

If no, make the following settings.

- Click on "New".
- Select "Completed".
- Select "Don't care".
- Select "Reports".
- Enter the AE title from PACS obtained from the customer under "to".
- Click on "OK".
- Close the windows.

**Note:** The required configuration for when the report is to come from PACS is now complete.



# DICOM Configuration of FLUOROSPOT Compact based systems (Luminos dRF, Luminos Agile, Ysio, UROSKOP Omnia)

CARE Analytics is freeware, a non-medical software tool designed to query dose information from Structured Dose Report objects stored on configured DICOM nodes. The tool can be executed from any PC in the hospital network or from a USB stick approved by a hospital IT administrator. CARE Analytics was tested with Structured Dose Reports from Luminos dRF, Luminos Agile, Ysio and UROSKOP Omnia VC10.

Strutured Dose Reports from FLUOROSPOT Compact based systems have to be send to a PACS system. CARE Analytics needs to be configured to retrieve the Structured Dose Reports from the PACS.

## Configuration of FLUOROSPOT Compact

### Step 1

Ask the customer for the following parameters.

Alias of DICOM store Node (PACS node)


### Step 2

**Configure the transfer of Structured Dose Reports to PACS node.**

Structured Dose Reports can be sent automatically to a PACS node, if a patient studie is finished, or sent manually when the images are sent. It is recommended to use either automatic or manual transfer.

In case of autotransfer of DICOM Dose SRs (shall be used, if patient studies are closed after examination):

- Select "User Settings"
- Select "DICOM Properties"
- Select "Sending".



- Enable "Send as DICOM Dose SR" in section "Autosending dose information" and select a destination (Alias from step 1)
- Click OK.


**Explanation:**

Now the reports are sent automatically, when the study is closed.



In case of manual transfer of DICOM Dose SRs (shall be used, if patient studies are not closed after examination)

- Select "User Settings"
- Select "DICOM Properties"
- Select "Sending".






- Enable "Send as DICOM Dose SR" in section "Consider exam protocol for manual send".
- Click Ok.

**Explanation:**

Now the reports could be send manually to the PACS together with the images. The checkbox DICOM Dose SR is enabled automatically.




# About CARE Analytics:

CARE Analytics is a freeware, non-medical software tool designed to extract dose information from Structured Dose Report objects stored on configured DICOM nodes or local media. The tool can be executed from any PC. The Structured Dose Report can be queried from any DICOM node within the hospital network or a USB stick approved by a hospital IT administrator. The Structured Dose Report objects are not made available on the PC, only distinct information from these objects is presented to the user or transformed into other formats. The information can be exported in xml-format to other programs like Microsoft® Office Excel for further analysis, e.g. statistical reports of dose data and related parameters.

CARE Analytics supports DICOM Structured Dose Reports from modality types:

- Computed Tomography
- Projection X-Ray
- Mammography

## General Layout



The screenshot displays the CARE Analytics software interface. The main window contains a large table with columns for PatientID, StudyDescription, SeriesDescription, SeriesTime, Procedure reported, System Name, Mean CTDvol, Total Number of Irradiation Events, Acquisition Protocol, Target Region, Scanning Length, Nominal Single Collimation Width, Pitch Factor, and Organ Characteristic. The table lists various CT scans for different patients, including head, spine, and abdomen scans. At the bottom of the window, there is a summary table with columns for DICOM Node, Study Date From, Study Date To, Patient Age From, Patient Age To, Patient ID, Performing Phys, Study Desc, System Name, Serial Number, Gender, Organ Characteristic, and Show Topogram Data. The summary table shows data for the year 2011-3-11, with a total of 16 patients, 16 studies, and 16 organs.

Gender	PatientID	StudyDescription	SeriesDescription	SeriesTime	Procedure reported	System Name	Mean CTDvol	Total Number of Irradiation Events	Acquisition Protocol	Target Region	Scanning Length	Nominal Single Collimation Width	Pitch Factor	Organ Characteristic
O	10.11.1...	Head*HeadRoutine (AdJ)	Dose Report	090329.32...	Computed Tomograph...	SOMATOM Def...	40.75 mGy	2 events	Head	Head	154 mm	0.6 mm	0.55 ratio	Head
O	10.11.1...	Head*HeadRoutine (AdJ)	Dose Report	153745.73...	Computed Tomograph...	SOMATOM Def...	36.69 mGy	2 events	Head	Head	154 mm	0.6 mm	0.55 ratio	Head
O	10.11.1...	Spine*SpineRoutine (AdJ)	Dose Report	154212.82...	Computed Tomograph...	SOMATOM Def...	8.17 mGy	2 events	Spine	Spine	366 mm	0.6 mm	0.8 ratio	Spine
O	10.11.1...	Head*HeadRoutine (AdJ)	Dose Report	154637.73...	Computed Tomograph...	SOMATOM Def...	266.03 mGy	2 events	Head	Head	143 mm	0.6 mm	0.35 ratio	Head
O	10.10.2...	Spine*SpineRoutine_Vol (AdJ)	Dose Report	120448.07...	Computed Tomograph...	SOMATOM Def...	2.36 mGy	3 events	Spine	Spine	300 mm	0.6 mm	0.8 ratio	Spine
O	10.10.2...	Thorax*ThoraxRt (AdJ)	Dose Report	144340.42...	Computed Tomograph...	SOMATOM Def...	3.36 mGy	2 events	Spine	Spine	200 mm	0.6 mm	0.8 ratio	Spine
O	10.10.2...	Spine*SpineRoutine_Vol (AdJ)	Dose Report	153537.42...	Computed Tomograph...	SOMATOM Def...	2.65 mGy	3 events	Spine	Spine	275 mm	0.6 mm	0.8 ratio	Spine
O	10.11.0...	Abdomen*AbdFullPhase (AdJ)	Dose Report	113510.75...	Computed Tomograph...	SOMATOM Def...	3.32 mGy	4 events	Non-Contrast	Abdomen	529 mm	1.2 mm	1.5 ratio	Abdomen
O	10.11.0...	Vascular*ChetECG_Abd (AdJ)	Dose Report	112141.78...	Computed Tomograph...	SOMATOM Def...	2.88 mGy	2 events	Arterial Phase	Abdomen	498 mm	0.6 mm	0.6 ratio	Abdomen
O	10.11.0...	Cardiac*FAST_Cardio (AdJ)	Dose Report	113617.75...	Computed Tomograph...	SOMATOM Def...	1.07 mGy	2 events	Venous Phase	Abdomen	498 mm	1.2 mm	0.6 ratio	Abdomen
O	10.11.0...	Cardiac*FAST_Cardio (AdJ)	Dose Report	113617.75...	Computed Tomograph...	SOMATOM Def...	0.48 mGy	5 events	Abdomen	Abdomen	282 mm	0.6 mm	0.6 ratio	Abdomen
O	10.11.0...	Cardiac*FAST_Cardio (AdJ)	Dose Report	113617.75...	Computed Tomograph...	SOMATOM Def...	17.79 mGy	5 events	CardSeq	Heart	104 mm	1.2 mm	0.9 ratio	Cardio
O	10.11.0...	Cardiac*FAST_Cardio (AdJ)	Dose Report	113617.75...	Computed Tomograph...	SOMATOM Def...	0.76 mGy	5 events	TestBkSeq	Chest	10 mm	10 mm	0.9 ratio	Thorax
O	10.11.0...	Cardiac*FAST_Cardio (AdJ)	Dose Report	113617.75...	Computed Tomograph...	SOMATOM Def...	1.21 mGy	5 events	CorCTA	Heart	104 mm	0.6 mm	0.9 ratio	Cardio
O	10.11.1...	Thorax*ThoraxRoutine (AdJ)	Dose Report	173940.81...	Computed Tomograph...	SOMATOM Def...	1.19 mGy	3 events	CorCTA	Heart	139 mm	0.6 mm	0.22 ratio	Cardio
O	10.11.1...	Head*HeadRoutine (AdJ)	Dose Report	174529.07...	Computed Tomograph...	SOMATOM Def...	40.75 mGy	2 events	PreflowMonitoring	Chest	10 mm	10 mm	1.2 ratio	Thorax
O	10.11.1...	Head*HeadRoutine (AdJ)	Dose Report	174529.07...	Computed Tomograph...	SOMATOM Def...	39.4 mGy	2 events	PreflowMonitoring	Chest	336 mm	0.6 mm	1.2 ratio	Thorax
O	10.11.1...	Special*PolyTrauma (AdJ)	Dose Report	175112.06...	Computed Tomograph...	SOMATOM Def...	8.35 mGy	5 events	Head	Head	154 mm	0.6 mm	0.55 ratio	Head
O	10.11.1...	Thorax*ThoraxRoutine (AdJ)	Dose Report	090358.50...	Computed Tomograph...	SOMATOM Def...	1.46 mGy	2 events	Neck	Neck	189 mm	1.2 mm	0.8 ratio	Neck
O	10.11.1...	Cardiac*FAST_Cardio (AdJ)	Dose Report	092104.32...	Computed Tomograph...	SOMATOM Def...	0.76 mGy	5 events	Neck	Neck	360 mm	1.2 mm	1.4 ratio	Abdomen
O	10.11.1...	Cardiac*FAST_Cardio (AdJ)	Dose Report	092104.32...	Computed Tomograph...	SOMATOM Def...	1.25 mGy	5 events	PreflowMonitoring	Chest	10 mm	10 mm	1.2 ratio	Thorax
M	10.03.1...	Shoulder*ShoulderRoutine (A...	Dose Report	091536.84...	Computed Tomograph...	SOMATOM Def...	1.18 mGy	2 events	CardSeq	Heart	104 mm	1.2 mm	0.9 ratio	Cardio
M	10.03.1...	Shoulder*ShoulderRoutine (A...	Dose Report	091536.84...	Computed Tomograph...	SOMATOM Def...	1.18 mGy	2 events	TestBkSeq	Chest	10 mm	10 mm	0.9 ratio	Thorax
M	10.03.1...	Shoulder*ShoulderRoutine (A...	Dose Report	091536.84...	Computed Tomograph...	SOMATOM Def...	1.18 mGy	2 events	CorAdSeq	Heart	104 mm	0.6 mm	0.9 ratio	Cardio
M	10.03.1...	Shoulder*ShoulderRoutine (A...	Dose Report	091536.84...	Computed Tomograph...	SOMATOM Def...	1.18 mGy	2 events	CorCTA	Heart	139 mm	0.6 mm	0.22 ratio	Cardio

DICOM Node	Study Date From	Study Date To	Patient Age From	Patient Age To	Patient ID	Performing Phys	Study Desc	System Name	Serial Number	Gender	Organ Characteristic	Show Topogram Data
	2011-3-11									ALL	N	

Figure 1: Main Screen of CARE Analytics (Computer Tomography)

## Workflow

- Configure a DICOM node in the <DICOM node configuration> dialog to query the Structured Dose Reports or import data directly from local media.
- Define Query Criteria and start the search or open an already saved file.
- Transfer the saved xml file to Microsoft® Office Excel for further evaluations.

# Toolbar

## 1. Query dose reports from a network node

Opens a window, where the user can define a query.

- Select the type of Structured Dose Report
- Enter the query criteria
- Select a preconfigured DICOM node
- Select "Search" Button to start the query from the defined DICOM node

### a. Dose SR type

Opens a drop down box, which allows the selection of SR types. It can be selected whether the query shall contain "Computed Tomography X-Ray", "Projection X-Ray", or "Mammography" dose reports.

### b. Query criteria

- Time Range <From> or <To> in YYYY-MM-DD
- Patient Birth Date <From> or <To> in years
- Patient ID (use wildcard character "\*" for broad search)
- Performing Physician
- Study Description
- System Name
- Serial Number of System
- Gender
- Organ Characteristic (only available for "Computed Tomography X-Ray" report)

Note: Characters compatible to ISO Latin-1 character sets can be used for query.

### c. Select a preconfigured DICOM node

Unless one or more DICOM nodes have been defined in the <DICOM node configuration> dialog, the drop-down box allows the selection of one predefined DICOM node for the query.

### d. Select "Search" Button to start the query from the defined DICOM node

The status bar will inform about the found reports and reports that were not considered for display due to the entered search criteria ("filtered out").




Figure 2: Toolbar




Figure 3: Query Criteria Dialog




Figure 4: Query Dose Reports from Local Media Dialog

## 2. Query dose reports from local media

Opens the <Query Criteria> dialog (as described above) containing an additional field "Import From..." to indicate the location of Structured Dose Report objects. The "Select" button brings up an "Open" dialog, which allows browsing through directories from local file systems.

The "Search" Button starts the import of the Structured Dose Report information into CARE Analytics according to the selected objects and entered criteria.

## 3. "Save file as xml"

A "Save as" dialog pops up where the filename and location can be entered. This file can be reopened for further evaluations.

## 4. Update dose reports

Starts a new query with the same query criteria. The new dose reports are added to the list.

## 5. Open file

Opens an already saved xml file. If the selected file is invalid, a message window pops up: "Invalid xml file". The status bar shows: "Unable to load report"

## 6. Merge files

Multiple existing XML files can be merged for further evaluations, e.g. in Microsoft® Office Excel

## 7. Export to Excel

In case Microsoft® Office Excel is installed on the same system as CARE Analytics, the displayed results can be exported to Microsoft® Office Excel for further analysis. Microsoft® Office Excel will open automatically when this function is selected.

The Microsoft® Office Excel output file includes different worksheets depending on the chosen report type "Computed Tomography" or "Projection X-Ray" or "Mammography":

### Computed Tomography

Select the desired template to be used for displaying the information in Microsoft® Office Excel.

Select the desired template to be used for displaying the information in Microsoft® Office Excel.

- Details  
The results are presented within Microsoft® Office Excel structured in two sheets, see also *Projection X-Ray / Mammography*
- Summary

Uses a template which includes 4 sheets:

### a. Selection Criteria

The query criteria are shown and the calculated max value for CTDI and DLP corresponding to every <Organ Characteristic>.

### b. CARE Analytics Data List

The Excel file contains the parameters listed in the CARE Analytics tool

### c. Ready-made Results

The Excel file contains dose related parameters for each Scan range:

- CTDIvol: {Min, Average, Max values}
- DLP: {Min, Average, Max values}
- mAs: {Min, Average, Max values}
- kV: {Min, Average, Max values}

### d. DLP Monitoring

The Excel file shows the DLP and CTDIvol distribution split per month for each organ characteristic

## Command Line Support

CARE Analytics comes with a command line version named "CATool.exe". The "CATool.exe" needs to be copied to the same folder as the files "CareAnalytics.exe" and "care\_conf.xml".

The command line support ensures CARE Analytics retrieves the required dose information automatically. Using the command line tool the user can easily define own queries without executing the interactive CARE

Analytics tool. It enables e.g. to write scripts for regular, automatic queries. The command line options are listed in figure 7 and can be printed on the screen using the "CATool.exe /h" command.

Example: "CATool.exe /sr ct /node NodeName" will retrieve CT dose report information from a preconfigured DICOM node called "NodeName".

```
D:\Programme\CareAnalytics\CARE_Analytics_U2_beta>catool.exe /?
Care Analytics Version U2
Copyright(C) Siemens AG 2011

Catool: query dose SR from DICOM node
usage: Catool.exe [options]

general options:
  /out File           File path of output xml result.
  /verbose            Print debug information.
  /h                 Print this help text and exit

query options:
  /sr SRTYPE          query dose SR type.
                     Possible type: CT, XRay, Mammo
  /node NodeName      Dicom node name which is configured in the GUI tool CareAnalytics.exe
  /sfrom YYYY-MM-DD   Lower bound of the range of study date.
  /sto YYYY-MM-DD     Upper bound of the range of study date.
  /bfrom YYYY-MM-DD   Lower bound of the range of patient's date of birth.
  /bto YYYY-MM-DD     Upper bound of the range of patient's date of birth.
  /pid PID            Patient ID.
  /pp PerformingPhysician Performing physician's name.
  /sd StudyDescriptio Description of the study.
  /sname SystemName   System Name. Only valid when SR type is CT.
  /sno SerialNumber   Serial Number.
  /g Gender            Possible value: M, F, O, A, MF
  /organ Organ        Organ characteristic. Only valid if the SR type is CT.
                     Possible value: Head, Neck, Shoulder, Thorax, Abdomen, Pelvis, Spine, Osteo, Extremitie
                     AngioHead, AngioBody, Runoff, PerfBody, Cardio, Respiratory
  /topo              Show tomogram info. Only valid if the SR type is CT.
  /y YYYY            Fetch all SR in the year YYYY. Only valid if either /sfrom or /sto is not used.
  /m MM              Fetch all SR from last MM months. Only valid if /sfrom, /sto or /y is not used.
  /w WW              Fetch all SR from last WW weeks. Only valid if /sfrom, /sto, /y or /m is not used.
```

Figure 8: Command Line Options

## Please mind

CARE Analytics is freeware. It might not be trained by Siemens application specialist and will not be installed by Siemens service engineers. The selection of a PC or local media within the hospital network is under the responsibility of the hospital IT administrator. The local

policies have to be observed to guarantee that there is no risk related to the attaching of this PC or USB stick to the DICOM network (i.e. virus protection). It is prohibited to copy CARE Analytics on any Siemens medical device.

## Toolbar

### *Computed Tomography Details / Projection X-Ray / Mammography*

Using Projection X-Ray or Mammography Structured Dose Reports the results are presented within Microsoft® Office Excel structured in two sheets:

#### **a. Patients**

This sheet contains all parameters selected in the CARE Analytics column configuration, which are either patient or cumulated dose information.

#### **b. Dose Events**

This sheet contains the information parameters selected in the CARE Analytics column configuration for each irradiation event.

### **8. Close current dose report**

Closes the current file. The user is asked whether he wants to save the file or not

### **9. Delete entry from the list**

Single or multiple entries (e.g. Phantom scans) can be deleted from the List

### **10. Column configuration**

The parameters which should be displayed into CARE Analytics and – in case of Computed Tomography Details view and “Projection X-Ray” reports – also transferred to Microsoft® Office Excel can be configured.

### **11. DICOM node configuration**

Multiple DICOM nodes can be configured. The default CARE Analytics configuration can be used for a single user.

### **12.Help**

Access to the “Help” file and “About CARE Analytics”.




Figure 5: Column Configuration Dialog (Computed Tomography)




Figure 6: DICOM Node Configuration




Figure 7: About Care Analytics Dialog



# COPYRIGHT File

This software is based in part on the work of the Independent JPEG Group. Also this software is partly based on CDragListBoxEx, which is licensed under The GNU Lesser General Public License (LGPLv3).

Unless otherwise specified, the DCMTK software package has the following copyright:

## Copyright (C) 1994-2004, OFFIS

This software and supporting documentation were developed by

Kuratorium OFFIS e.V.  
Healthcare Information and Communication Systems  
Escherweg 2  
D-26121 Oldenburg, Germany

THIS SOFTWARE IS MADE AVAILABLE, AS IS, AND OFFIS MAKES NO WARRANTY REGARDING THE SOFTWARE, ITS PERFORMANCE, ITS MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR USE, FREEDOM FROM ANY COMPUTER DISEASES OR ITS CONFORMITY TO ANY SPECIFICATION. THE ENTIRE RISK AS TO QUALITY AND PERFORMANCE OF THE SOFTWARE IS WITH THE USER.

Copyright of the software and supporting documentation is, unless otherwise stated, owned by OFFIS, and free access is hereby granted as a license to use this software, copy this software and prepare derivative works based upon this software. However, any distribution of this software source code or supporting documentation or derivative works (source code and supporting documentation) must include the three paragraphs of this copyright notice.

Some portions of the DCMTK software package are derived from earlier versions of this software with the following copyright, and can be identified by the following copyright notice located in each source file:

## Copyright (C) 1993/1994, OFFIS, Oldenburg University and CERIUUM

This software and supporting documentation were developed by

Institut OFFIS  
Bereich Kommunikationssysteme  
Westerstr. 10-12  
26121 Oldenburg, Germany

Fachbereich Informatik  
Abteilung Prozessinformatik  
Carl von Ossietzky Universitaet  
Oldenburg  
Ammerlaender Heerstr. 114-118  
26111 Oldenburg, Germany

CERIUUM  
Laboratoire SIM  
Faculte de Medecine  
2 Avenue du Pr. Leon Bernard  
35043 Rennes Cedex, France

for CEN/TC251/WG4 as a contribution to the Radiological Society of North America (RSNA) 1993 Digital Imaging and Communications in Medicine (DICOM) Demonstration.

THIS SOFTWARE IS MADE AVAILABLE, AS IS, AND NEITHER OFFIS, OLDENBURG UNIVERSITY NOR CERIUUM MAKE ANY WARRANTY REGARDING THE SOFTWARE, ITS PERFORMANCE, ITS MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR USE, FREEDOM FROM ANY COMPUTER DISEASES OR ITS CONFORMITY TO ANY SPECIFICATION. THE ENTIRE RISK AS TO QUALITY AND PERFORMANCE OF THE SOFTWARE IS WITH THE USER.

Copyright of the software and supporting documentation is, unless otherwise stated, jointly owned by OFFIS, Oldenburg University and CERIUUM and free access is hereby granted as a license to use this software, copy this software and prepare derivative works based upon this software. However, any distribution of this software source code or supporting documentation or derivative works (source code and supporting documentation) must include the three paragraphs of this copyright notice.

Some other parts of this software within the dcmstk/dcmnet sub-package related to the DICOM Upper Layer Protocol are derived from software developed for the RSNA'93 DICOM demonstration and kindly made available to us by the Mallinckrodt Institute of Radiology. Such software can be identified by the following copyright notice located in each affected source file:

## Copyright (C) 1993, RSNA and Washington University

The software and supporting documentation for the Radiological Society of North America (RSNA) 1993 Digital Imaging and Communications in Medicine (DICOM) Demonstration were developed at the

Electronic Radiology Laboratory  
Mallinckrodt Institute of Radiology  
Washington University School of Medicine  
510 S. Kingshighway Blvd.  
St. Louis, MO 63110, USA

as part of the 1993 DICOM Central Test Node project for, and under contract with, the Radiological Society of North America.

THIS SOFTWARE IS MADE AVAILABLE, AS IS, AND NEITHER RSNA NOR WASHINGTON UNIVERSITY MAKE ANY WARRANTY ABOUT THE SOFTWARE, ITS PERFORMANCE, ITS MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR USE, FREEDOM FROM ANY COMPUTER DISEASES OR ITS CONFORMITY TO ANY SPECIFICATION. THE ENTIRE RISK AS TO QUALITY AND PERFORMANCE OF THE SOFTWARE IS WITH THE USER.

Copyright of the software and supporting documentation is jointly owned by RSNA and Washington University, and free access is hereby granted as a license to use this software, copy this software and prepare derivative works based upon this software. However, any distribution of this software source code or supporting documentation or derivative works (source code and supporting documentation) must include the three paragraphs of the copyright notice.

The dcmjpeg sub-package includes an adapted version of the Independent JPEG Group Toolkit Version 6b, which is contained in dcmjpeg/libijg8, dcmjpeg/libijg12 and dcmjpeg/libijg16. This toolkit is covered by the following copyright. The original README file for the Independent JPEG Group Toolkit is located in dcmjpeg/docs/ijg\_readme.txt.

## dcmjpeg Sub-Package

The authors make NO WARRANTY or representation, either express or implied, with respect to this software, its quality, accuracy, merchantability, or fitness for a particular purpose. This software is provided "AS IS", and you, its user, assume the entire risk as to its quality and accuracy.

This software is copyright (C) 1991-1998, Thomas G. Lane. All Rights Reserved except as specified below.

Permission is hereby granted to use, copy, modify, and distribute this software (or portions thereof) for any purpose, without fee, subject to these conditions:

1. If any part of the source code for this software is distributed, then this README file must be included, with this copyright and no-warranty notice unaltered; and any additions, deletions, or changes to the original files must be clearly indicated in accompanying documentation.
2. If only executable code is distributed, then the accompanying documentation must state that "this software is based in part on the work of the Independent JPEG Group".
3. Permission for use of this software is granted only if the user accepts full responsibility for any undesirable consequences; the authors accept NO LIABILITY for damages of any kind.

These conditions apply to any software derived from or based on the IJG code, not just to the unmodified library. If you use our work, you ought to acknowledge us.

Permission is NOT granted for the use of any IJG author's name or company name in advertising or publicity relating to this software or products derived from it. This software may be referred to only as "the Independent JPEG Group's software".

We specifically permit and encourage the use of this software as the basis of commercial products, provided that all warranty or liability claims are assumed by the product vendor.

The color quantization code in module `dcmimage` (`dcmquant` and the related classes) is derived from code written by Jef Poskanzer for the NetPBM toolkit which has the following copyright:

Copyright (C) 1989, 1991 by Jef Poskanzer.

Permission to use, copy, modify, and distribute this software and its documentation for any purpose and without fee is hereby granted, provided that the above copyright notice appear in all copies and that both that copyright notice and this permission notice appear in supporting documentation. This software is provided "as is" without express or implied warranty.

The code for the `OFStandard::strlcpy` and `OFStandard::strlcat` helper functions in `ofstd/libsrc/ofstd.cc` has been derived from the BSD implementation of `strlcpy()` and `strlcat()` and which carries the following copyright notice:

Copyright (c) 1998 Todd C. Miller <Todd.Miller(at)courtesan.com>  
All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

1. Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
2. Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
3. The name of the author may not be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE AUTHOR BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

The code for the `OFStandard::atof` helper function in `ofstd/libsrc/ofstd.cc` has been derived from an implementation which carries the following copyright notice:

## Copyright 1988 Regents of the University of California

Permission to use, copy, modify, and distribute this software and its documentation for any purpose and without fee is hereby granted, provided that the above copyright notice appear in all copies. The University of California makes no representations about the suitability of this software for any purpose. It is provided „as is“ without express or implied warranty.

The code for `OFStandard::ftoa` has been derived from an implementation which carries the following copyright notice:

## Copyright (c) 1988 Regents of the University of California. All rights reserved.

Redistribution and use in source and binary forms are permitted provided that the above copyright notice and this paragraph are duplicated in all such forms and that any documentation, advertising materials, and other materials related to such distribution and use acknowledge that the software was developed by the University of California, Berkeley. The name of the University may not be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED "AS IS" AND WITHOUT ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, WITHOUT LIMITATION, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

The "Base64" encoder/decoder in `ofstd/libsrc/ofstd.cc` has been derived from an implementation which carries the following copyright notice:

## Copyright (c) 1999, Bob Withers – [bwit\(at\)pobox.com](mailto:bwit(at)pobox.com)

This code may be freely used for any purpose, either personal or commercial, provided the authors copyright notice remains intact.

The dcmjp2k sub-package (which is currently not part of the free toolkit) includes an adapted version of the JasPer JPEG 2000 toolkit, which is contained in dcmjp2k/libjaspr. This toolkit is covered by the following copyright.

## JasPer License Version 2.0

Copyright (c) 1999-2000 Image Power, Inc.  
Copyright (c) 1999-2000 The University of British Columbia  
Copyright (c) 2001-2003 Michael David Adams

All rights reserved.

Permission is hereby granted, free of charge, to any person (the "User") obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

1. The above copyright notices and this permission notice (which includes the disclaimer below) shall be included in all copies or substantial portions of the Software.
2. The name of a copyright holder shall not be used to endorse or promote products derived from the Software without specific prior written permission.

THIS DISCLAIMER OF WARRANTY CONSTITUTES AN ESSENTIAL PART OF THIS LICENSE. NO USE OF THE SOFTWARE IS AUTHORIZED HEREUNDER EXCEPT UNDER THIS DISCLAIMER. THE SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT OF THIRD PARTY RIGHTS. IN NO EVENT SHALL THE COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, OR ANY SPECIAL INDIRECT OR CONSEQUENTIAL DAMAGES, OR ANY DAMAGES WHATSOEVER RESULTING FROM LOSS OF USE, DATA OR PROFITS, WHETHER IN AN ACTION OF CONTRACT, NEGLIGENCE OR OTHER TORTIOUS ACTION, ARISING OUT OF OR IN CONNECTION WITH THE USE OR PERFORMANCE OF THIS SOFTWARE. NO ASSURANCES ARE PROVIDED BY THE COPYRIGHT HOLDERS THAT THE SOFTWARE DOES NOT INFRINGE THE PATENT OR OTHER INTELLECTUAL PROPERTY RIGHTS OF ANY OTHER ENTITY. EACH COPYRIGHT HOLDER DISCLAIMS ANY LIABILITY TO THE USER FOR CLAIMS BROUGHT BY ANY OTHER ENTITY BASED ON INFRINGEMENT OF INTELLECTUAL PROPERTY RIGHTS OR OTHERWISE. AS A CONDITION TO EXERCISING THE RIGHTS GRANTED HEREUNDER, EACH USER HEREBY ASSUMES SOLE RESPONSIBILITY TO SECURE ANY OTHER INTELLECTUAL PROPERTY RIGHTS NEEDED, IF ANY. THE SOFTWARE IS NOT FAULT-TOLERANT AND IS NOT INTENDED FOR USE IN MISSION-CRITICAL SYSTEMS, SUCH AS THOSE USED IN THE OPERATION OF NUCLEAR FACILITIES, AIRCRAFT NAVIGATION OR COMMUNICATION SYSTEMS, AIR TRAFFIC CONTROL SYSTEMS, DIRECT LIFE SUPPORT MACHINES, OR WEAPONS SYSTEMS, IN WHICH THE FAILURE OF THE SOFTWARE OR SYSTEM COULD LEAD DIRECTLY TO DEATH, PERSONAL INJURY, OR SEVERE PHYSICAL OR ENVIRONMENTAL DAMAGE ("HIGH RISK ACTIVITIES"). THE COPYRIGHT HOLDERS SPECIFICALLY DISCLAIM ANY EXPRESS OR IMPLIED WARRANTY OF FITNESS FOR HIGH RISK ACTIVITIES.

# DICOM Information

CARE Analytics Tool provides standard conformance to the following DICOM 3.0 SOP Classes as a SCU:

Verification SOP Class	1.2.840.10008.1.1
Study Root Query/Retrieve Info Model – FIND	1.2.840.10008.5.1.4.1.2.2.1
Study Root Query/Retrieve Info Model – MOVE	1.2.840.10008.5.1.4.1.2.2.2

The transfer syntax is DICOM Implicit VR Little Endian (1.2.840.10008.1.2)

CARE Analytics Tool provides standard conformance to the following DICOM 3.0 SOP Classes as a SCP  
(Default AE title: CareAnalytics, Default IP port: 5104)

Verification SOP Class	1.2.840.10008.1.1
X-Ray Radiation Dose Structured Reports Storage	1.2.840.10008.5.1.4.1.1.88.67

The transfer syntax is DICOM Implicit VR Little Endian (1.2.840.10008.1.2)



#### **Global Business Unit**

Siemens AG  
Medical Solutions  
Computed Tomography & Radiation Oncology  
Siemensstr. 1  
DE-91301 Forchheim  
Germany  
Phone: +49 9191 18-0  
Fax: +49 9191 18 9998

Siemens AG  
Medical Solutions  
Angiography & Interventional X-Ray Systems  
Siemensstr. 1  
DE-91301 Forchheim  
Germany  
Phone: +49 9191 18-0  
Fax: +49 9191 18 9998

Siemens AG  
Medical Solutions  
X-Ray Products  
Henkestr. 127  
DE-91052 Erlangen  
Germany  
Phone: +49 9131 84-0

#### **Global Siemens Headquarters**

Siemens AG  
Wittelsbacherplatz 2  
80333 Muenchen  
Germany

#### **Global Siemens Healthcare Headquarters**

Siemens AG  
Healthcare Sector  
Henkestr. 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0  
[www.siemens.com/healthcare](http://www.siemens.com/healthcare)

#### **Legal Manufacturer**

Siemens AG  
Wittelsbacherplatz 2  
DE-80333 Muenchen  
Germany

[www.siemens.com/healthcare](http://www.siemens.com/healthcare)